

The correct handling and use of the Wincorn Medical Aluminum Hand is described in the following.


In order to keep risks and unnecessary burdens as low as possible for the patient, the user, and third parties, we request that you carefully look through these instructions for use and keep them for further reference.


1. Intended use

The Wincorn Medical aluminum hand can be used in all surgical interventions in the hand region. It is intended to fixate and immobilize the hand to be operated on. The hand to be operated on, and the patient's individual fingers, are fixated and immobilized by bending the wrist tabs and fingertips of the aluminum hand over the patient's wrist and individual fingers.

- Place the patient's hand on top of the surgical aluminum hand.
- Bend wrist tabs over the patient's wrist and fingertips over the patient's fingers as needed. The patient's hand should then be flat and immobilized.
- The tab opposite of the thumb can be used as a "kickstand" to prevent the wrist from rolling.
- After the procedure, straighten the fingertips and wrist tabs on the aluminum hand.

2. Care and Cleaning

 The product is delivered non-sterile. The procedures outlined below should be followed to ensure the safe handling of biologically contaminated aluminum hands. **Sterilize prior to use!**

 Prior to each use, inspect the aluminum hand for any broken, worn, cracked, or fractured components. If any components are damaged or **if the black nylon coating is cracked**, discard the aluminum hand! **DO NOT USE!** Handle and dispose of in accordance with accepted medical practice and applicable laws and regulations.

The aluminum hand must be cleaned prior to sterilization in accordance with standard institutional protocol. Any non-acidic, general-purpose cleaner suited for standard surgical instruments is acceptable.

After cleaning, the aluminum hand should then be autoclaved. The aluminum hand is a reusable device designed to be used with the following sterilization cycle parameters:

Pre Vacuum cycle:

Double-wrapped in standard central supply wrap

for 4 minutes

132°C (270°F)

Drying Time minimum 4 minutes

Gravity cycle:

Double-wrapped in standard central supply wrap

for 15-30 minutes

121°C (250°F)

Drying Time minimum 15-30 minutes

Store in a dry, clean, and safe place.

Wincorn Medical does not accept any liability if these instructions for use have not been observed or followed.

Replacement parts can be obtained from Wincorn Medical:
11.2014 Aluminum Hand Adult 14" with tab, nylon coating, black

Meaning of the Symbols



Attention, observe notes.



The product is delivered non-sterile.



Reference number of an item

LOT Identification number assigned to a particular quantity or lot of material from the manufacturer.



Read the instructions for use.



Manufacturer's declaration that products comply with requirements of the European Commission Product Directives.



Manufacturer:

Address: 100 N Central Expy Suite 1310 Chase Bank Building Richardson, TX 75080

Phone: +1 214-630-1221

Email: info@kgwlawfirm.com

Website: www.kgwlawfirm.com